


TempPro


Fluorescence lifetime system

Fluorescence lifetime spectroscopy for every laboratory

TempPro brings the power of lifetime spectrofluorometry to researchers who previously found it unaffordable. Recording fluorescence lifetimes gives you a dynamic picture of your sample's characteristics, including molecular motion and size, solvation environment, and intermolecular distances—things that steady-state fluorescence cannot provide. Lifetime spectrofluorometry is unique in its sensitivity and selectivity: for example, you can temporally separate similar fluorescence spectra from two distinct species in a solution. TempPro uses Time-Correlated Single-Photon Counting (TCSPC)—the reliable, intuitive, and most sensitive lifetime-measurement technique. The filter-based TempPro can measure a complete range of lifetimes from only 100 picoseconds to 1 second, with our NanoLED and SpectraLED solid-state pulsed excitation sources spanning wavelengths from the UV to NIR. The TempPro lifetime spectrofluorometer is so affordable that any lab now can exploit the power of fluorescence dynamics coupled with a phosphorescence capability.



Fluorescence decay of 10^{-5} M aqueous Erythrosin B, taken with the TempPro, giving a lifetime of 85 ps. Excitation = 488 nm. Emission > 550 nm. The green line is the mono-exponential fit to the decay data (red dots) after deconvolution with the instrumental response (blue dots).



Anisotropy decay of 10^{-5} M DPH in Kaydol®, recorded with the TempPro. Excitation = 373 nm. Emission > 418 nm. The initial anisotropy of 0.33 drops to zero showing a free rotation with 4.9 ns correlation time.

APPLICATIONS


- FRET (Förster Resonance Energy Transfer)
- Stern-Volmer quenching
- Efficiency of photosynthesis
- Lanthanide luminescence
- Fluorescence anisotropy
- Protein fluorescence
- Phosphorescence lifetimes and anisotropy



Full performance at a desktop price!

Value, power, and compactness are the hallmarks of TemPro, a product perfect for basic research and analytical measurements. HORIBA Scientific backs TemPro with nearly 200 years of sales, service, and applications expertise in optical instrumentation.


ADVANCED FEATURES



Span 100 ps to 1 s with the twist of a dial using TemPro's intuitive data-acquisition and lifetime-analysis software, DataStation and DAS6.

- Lifetime determination from 100 ps to 1 s*
- Detector response from 185–650 nm (others optional)
- Single-photon counting detection
- Data-acquisition and analysis software
- Fast, reliable USB 2.0 data link

* Dependent on source, sample, and detector.



Upgrades

- NIR extended range detector (to 850 nm)
- Digital temperature sensor and magnetic stirring
- Front-face sample holder
- Monochromators
- Polarizers

Our versatile pulsed solid-state excitation sources: NanoLED and SpectraLED.



HORIBA Scientific provides reliability and expertise at a surprising price!

info-sci@horiba.com
www.horiba.com/scientific

HORIBA
Scientific

USA: +1 732 494 8660
UK: +44 (0)20 8204 8142
Spain: +34 91 490 23 34
Other Countries: +33 (0)1 64 54 13 00
France: +33 (0)1 64 54 13 00
Italy: +39 0 2 5760 3050
China: +86 (0)10 8567 9966

Germany: +49 (0)89 4623 17-0
Japan: +81 (0)3 38618231
Brazil: +55 11 5545 1540

