

Measurement of Sulfur in Rubber

Alain Salaville, Jérôme Barraqué, Application Laboratory, Longjumeau, FRANCE

Keywords: rubber, sulfur, EMIA 820V

1 Introduction

Sulphur is added in the natural rubber to produce plastic materials. However, it is important to measure the sulphur contained in the rubber to define the chemical and mechanical properties of the final product.

concentration.

2 Instrumentation

2.1 Principle

The test was performed on the model EMIA 820V. The measurement principle is shown in Figure 2.

The sample is placed in a ceramic crucible in a

2.2 Unique Features

2.2.1 - Programmable Temperature Curves

high frequency induction furnace. The sample is heated at a programmable temperature. Gases produced during the combustion are then ana-

lyzed using four Infrared detectors, after dust

and moisture removal. The analysis of SO₂

determines sulfur concentration. The analysis of

low and high CO₂ and CO determine carbon

The high frequency or induction furnace is equipped with a plate current control function. This allows users to easily optimize the temperature according to the samples. Some customized temperature curves can be created in order to observe various phenomena such as surface con-











tamination and different phases or forms of carbon and sulfur.

2.2.2 - Direct gas analysis without conversion Four Infrared analyzers (NDIR) are used to directly analyze CO, CO₂ and SO₂ over the full range of concentrations. No converter is used nor cellulose filter to trap SO₃ generated in the converter.

2.2.3 - Computer System

All EMIA Series Analyzers are operated by a separate computer system. The software is compatible with Windows 95/98/2000/NT/XP. It includes several functions such as maintenance, diagnostis, statistical studies, curve and data traceability, etc.

2.2.4 - Automatic Cleaning

The double Auto Cleaner option features two brushes to simultaneously clean the combustion tube and the cylindrical dust filter after each measurement. The dust is removed to the dust box by a difference in pressure, which avoids the need for an external vacuum cleaner.

2.2.5 - Automation

It is possible to add standard modules for partial to full automation for 24/7 operation. For more detail see EA.TN 26: Options for Partial and Complete Automation.




Figure 2: Operating principle

3 Sample preparation

The sample was in the form of a powder.

- 1. Weight 0.05 g of sample into a ceramics crucible
- 2. Weight 0.5 g of pure iron, 1.5g of Tungsten and 0.3 g Tin as accelerator, And cover the sample with each metal.

4. Set the ceramic crucible with sample on the crucible stand, and press the [START] button to start analysis.



4 Conditions of analysis

Table 1: Operating conditions

p	Start oower mA)	End power (mA)	Time from start to end power (sec)
Step 1	0	175	5
Step 2 1	75	175	30

	Carbon	Sulfur
Purge time		15 sec
Integration wait time		5 sec
Integration time		100 sec
Comparator level		2.0 %
Comparator wait time	:	40 sec

5 Calibration

- 1. Set up the system to the analytical condition for the steel in the operator's instruction manual.
- 2. Calibrate the system following the procedure in the operator's instruction manual.
- 3. Weight 1.5g of Tungsten and 0.3g of Tin as blank into a ceramics crucible. Enter 1.0g as sample weight for blank analysis. Repeat measurement 3 times at minimum.
- 4. Weight 1.0g of JSS 243-4 (S:0.346mass%) into a ceramics crucible. And cover the sample with 1.5g of Tungsten and 0.3g of Tin. Repeat measurement 3 times at minimum.
- 5. Change sample analysis condition to the above table condition.
- 6. Compensate the blank signal because analytical condition to steel standard sample and Rubber sample is different. (As for the details, refer to the content of the blank shift of the instruction manual.)

7. Weight 0.5g of Pure Iron, 1.5g of Tungsten and 0.3g of Tin in the crucible. Enter 0.05g as sample weight for blank analysis. Repeat measurement 3 times at minimum.

6 Results on rubber carbide

Table 2: Rubber

Weight (g) Carbon (mass%)) Sulfur (mass%)
0.05102	3.159
0.05052	3.122
0.05087	3.127
0.05112	3.174
0.05092	3.170
Average	3.150
Standard Deviation	0.024
RSD(%)	0.77
Range	0.052



7 Summary

Instrument: EMIA-820V C/S Determinator

Calibration: JSS 243-4 (S: 0.346 mass %) 1.0 g

Sample: Rubber

Type: Chip Weight: 0.05 g

Accelerator: Pure iron (P/N 905.110.300.001) 0.5g

Tungsten (P/N 905.110.140.001) 1.5g Tin (P/N 905.202.200.001)0.3g

Crucible: Ceramic (P/N 905.202.200.001)

8 Conclusion

Sulfur measurement in rubber samples is compatible with the EMIA 820 V Series equipped with a high frequency furnace. The extraction is complete and efficient in all cases, and the results are repeatable.

This technical note is adapted from an Horiba technical note.

In the USA: HORIBA Jobin Yvon Inc. 3880 Park Avenue Edison, NJ 08820 Tel: 1-732-494-8660 Fax: 1-732-494-8796 E-mail: info@jobinyvon.com In France: HORIBA Jobin Yvon S.A.S. 16-18, rue du Canal 91165 Longjumeau Cedex Tel: (33) 1/64 54 13 00 Fax: (33) 1/69 09 90 88

Germany: (49) 89/46 23 17-0 Italy: (39) 2/57 60 56 90 U.K.: (44) 20/82 04 81 42 In Japan: HORIBA Ltd. 2 Miyanohigashi, Kisshoin Minami-ku, Kyoto 601-8510 TEL: (81) 75 313 8121 FAX: (81) 75 321 5725 www.jyhoriba.jp

China: (86) 10/6849 22 16 Spain: (34) 91/724 16 57 Other Countries: Contact HJY S.A.S.

HORIBAJOBIN YVON

www.jobinyvon.com