


Ultima Expert LT

High performance at affordable price



Designed for Challenging Applications

HORIBA Scientific develops and manufactures high performance ICP-OES spectrometers for more than 35 years.

The new Ultima Expert LT provides high performance at affordable price for laboratories with challenging samples

Ultima Expert LT integrates high efficiency Jobin Yvon optical design capable to achieve optimal performance for a large variety of sample types and matrices.

Ultima Expert LT is driven by the powerful Analyst software featuring a large variety of analytical functionalities for tailored control and analysis.

The robustness of the Ultima Expert LT makes it ideal for applications common to mining, chemicals manufacture, salt production, wear metals in oil analysis, petrochemical and metallurgical production

Gain in performance with the Ultima Expert LT for your most challenging applications!

Ultima Expert LT

Superior Performance in ICP-OES

Unique Plasma Torch Design for Most Versatile and Accurate Analyses


Ultima Expert LT uses a unique plasma torch design with radial viewing mode to allow viewing of the entire normal analytical zone.

The vertical torch, the original sheath gas device and the wide injector enable the Ultima Expert LT to handle difficult matrices, with robust operation and minimal maintenance.

Ultima Expert LT ensures the lowest detection limits even with the most challenging samples, such as high salt content, brines, dissolved solids and complex organics.

Its water-cooled high efficiency 40.68 MHz solid state generator enhances the ability to handle a large range of samples, reduces warm-up time down to 15 minutes, and improves stability and reliability.

All these features make the Ultima Expert LT the most versatile, accurate and robust ICP OES spectrometer available.



High Resolution and Full Wavelength Coverage

Ultima Expert LT delivers the highest resolution with less than 5 picometers for the UV range and less than 10 picometers for the Visible range. Such resolution can be achieved due to the unique optical design of the Ultima Expert LT, which integrates a high density holographic grating and one meter focal length optics.

Full wavelength coverage from 120 to 800 nm is offered to satisfy all the requirements for elemental analysis, including far UV capability for halogen elements analysis.

High resolution and excellent sensitivity allow elements with high, low and trace concentrations to be measured accurately, giving you maximum confidence in your results.

Easy to Handle, Easy to Maintain

Spacious compartments for sample and plasma are designed to facilitate sample handling.


Quick release, fully demountable torch with no adjustment required facilitates multi-user operation and provides excellent reproducibility.

Powerful Software with advanced features

The Analyst software of the Ultima Expert LT extends instrument capability beyond analytical performance with advanced feature for the user.

Complemented with S³-base viewer for the display of the information contained in the S³-base database for facilitated identification of potential spectral interferences, and an uncertainty calculation tool for calibration quality improvement, the Analyst software provides full instrument control and advanced features for method development, analysis and results management.

Thanks to all the functions available in Analyst software, Ultima Expert LT delivers results you can trust!



Ultima Expert LT

Specifications



Standard Configuration

Generator	Radio-frequency, solid-state 40.68 MHz, water-cooled
Spectral range	120 – 800 nm
Spectrometer	Thermally stabilized, 1 meter focal length with 2400 g/mm grating used in the 1st and 2nd order < 5 µm for 120-320 nm and < 10 µm for 320-800 nm
Resolution	Fully demountable torch with 3 mm i.d. alumina injector and quartz tubes
Plasma torch	Concentric glass nebulizer and glass cyclonic spray chamber, 3 channel peristaltic pump
Sample introduction	

Options

Instrument	Far UV kit for sensitive measurement for halogen elements
Accessories	Autosampler AS-500 with optional rinse station Argon humidifier Introduction system kits for improved performance (small volume, organics, Hydrofluoric acid, high total dissolved solids) Concomitant Metals Analyzer for simultaneous measurement of hydride forming elements and other elements Oxygen kit for alkali elements in organics

Facility requirements

Dimension (wxdxh)	1696 x 698 x 604 mm
Weight	205 kg (452 lb)
Power supply	Single phase, 220-240 V, 50-60 Hz, 4 kVA
Environmental	20 to 80 % humidity, 18-24°C at ± 2°C
Argon	99.995 % purity
Nitrogen	160 to 190 nm, 99.999 % purity 120 to 160 nm, 99.9995 % purity
Exhaust	250 m³/h (150 cfm)

www.horiba.com/scientific
info.sci@horiba.com

France:	HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel: +33 (0)1 69 74 72 00 - Fax: +33 (0)1 69 09 07 21 - Email: info-sci.fr@horiba.com
USA:	HORIBA Instruments Inc., 3880 Park Avenue, Edison, NJ 08820-3012 - Toll-free: +1-866-562-4698 - Tel: +1 732 494 8660 - Fax: +1 732 549 5125 - Email: info-sci.us@horiba.com
Japan:	HORIBA Ltd., Tokyo Branch Office, 2-6, KandaAwaji-cho, Chiyoda-ku, Tokyo 101-0063, Japan - Tel: +81-(0)3 6206 4721 - Fax: +81 (0)3 6206 4730 - Email: info-sci.jp@horiba.com
Germany:	HORIBA Jobin Yvon GmbH, Hauptstrasse 1, 82008 Unterhaching - Tel: +49 (0)89 4623 17-0 - Fax: +49 (0)89 4623 17-99 - Email: info-sci.de@horiba.com
Italy:	HORIBA Jobin Yvon Srl., Via Cesare Pavese 21, 20090 Opera (Milano) - Tel: +39 2 5760 3050 - Fax: +39 2 5760 0876 - Email: info-sci.it@horiba.com
UK:	HORIBA UK Ltd., 2 Dalston Gardens, Stanmore, Middlesex HA7 1BQ - Tel: +44 (0)20 8204 8142 - Fax: +44 (0)20 8204 6142 - Email: info-sci.uk@horiba.com
China:	HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel: +86 (0)21 6289 6060 - Fax: +86 (0)21 6289 5553 Email: info-sci.cn@horiba.com
Brazil:	HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multivias, Bairro Medeiros, Jundiaí / SP, CEP 13.212-181 - Tel: +55 (0)11 2923 5400 Fax: +55 (0)11 2923 5490 - Email: infoconfitecnica.br@horiba.com
Other:	Tel: +33 (0)1 69 74 72 00 - Email: info.sci@horiba.com