

**Your new business partner for
RoHS, Halogen free, ELV applications**



MESA-50



X-RAY FLUORESCENCE ANALYZER



Touch and

HORIBA has been selling our XGT-WR series of EDXRF analyzers worldwide for many years to provide screening measurements of samples containing hazardous elements such as Pb, Cd, Cr, Hg, Br for RoHS, ELV, and Cl for halogen-free applications.

XGT systems are used daily throughout the world.

HORIBA has now developed the new MESA-50 EDXRF analyzer, based on our vast experience with customer requirements and knowledge. MESA-50 provides user friendly operation and excellent performance. MESA-50 includes three analysis diameters, suitable for every sample, from thin cables and electronic parts to bulk samples. The unique combination of SDD (Silicon Drift Detector) and DPP (Digital Pulse Processor) changes the landscape of EDXRF analysis.

HORIBA's new MESA-50 supports ecological procurement; it contributes not only to EU RoHS and ELV compliance testing, but also regulatory work for many other countries.



Outstanding 5S Technology



1. Speedy

- **SDD (Silicon Drift Detector) drastically reduces measurement times, and provides increased sensitivity, for true high throughput analysis.**

2. Small

- **Portable, small footprint and light weight**
- **Internal battery power supply**

3. Simple

- **Reduces routine maintenance work (LN₂ free operation)**
- **No need for vacuum pumps**
- **Intuitive simple measurement process for all material types**

4. Smart

- **English/Japanese/Chinese user interfaces**
- **Excel® data management tool**

5. Safe


- **No worry about X-ray leakage**

Feel!

MESA-50

Hg Cr Cd Pb Br Cl Sb As

MESA-50 optics




Blue: MESA-50 $\phi 1.2$ mm, Red: Previous model (XGT-1000WR $\phi 1.2$ mm)

MESA-50 with SDD can:

- 1) Increase spectrum intensity and resolution
- 2) Decrease background

Compact optical design makes MESA-50 a portable instrument. The small X-ray generator with low power consumption permits the battery powered mode.


Better performance with special filters



Blue: MESA-50 $\phi 1.2$ mm

The X-ray tube target material is optimized for high sensitivity Cl analysis.

The MESA-50 doesn't require a vacuum pump for Cl measurement.



Yellow: MESA-50 $\phi 3$ mm, Red: Previous model (XGT-1700WR $\phi 3$ mm)

3 size collimators

$\phi 1.2$ mm, 3 mm, 7 mm (Automatic switching)


Operator can select the most suitable collimator size. The precise analysis area is displayed on the optical image. The small 1.2 mm beam means that operators can simply analyze individual cables, without needing to line up multiple cables to cover the X-ray beam area.



Simple operation, excellent support tools

1. Position the sample


The clear optical image from the CCD camera allows easy positioning of the sample. The visualized measurement point is shown on the CCD image.



2. Click the measurement button

Optimized measurement conditions for each material type are assigned to each measurement button (maximum 18 buttons). One click starts the measurement.


Operator does not need to set conditions each time, and no filter switching is necessary, thus minimizing possible operator error. Customized measurement buttons can be created, with original settings.



3. Display of spectrum and calculation result

Automatic display of quantitative results (concentration & 3 sigma) after measurement.

Spectrum




Result

Element	Conc.	3σ	Unit
Cr	1.00	0.05	ppm
Zn	1.00	0.05	ppm
Pb	1.00	0.05	ppm
Sr	1.00	0.05	ppm
Cd	1.00	0.05	ppm
Sn	1.00	0.05	ppm

4. Excel® data management software and report format

Automatically transfer to Excel® with one button click. If internal administration values are registered, OK/NG (No Good) judgments (results vs threshold values) are performed. Measurement report is also made with one button click.



Seq.	Part #	Part Name	Operator	Date	Part	Cr	Pb	Cd	Sn
1	0A0001	Cable1	HORIBA	07/20/12	Center	13.9	17.7	OK	125.2
2	0A0002	Cable2	HORIBA	07/20/12	Center	0.0	8.0	OK	112.5
3	0A0003	Cable3	HORIBA	07/20/12	Center	3.0	3.3	OK	132.5
4	0A0004	Cable4	HORIBA	07/20/12	Center	1.9	2.0	ND	125.1

Portable! Small size and battery system



Less than One Cubic Foot!

Size: 208(W) x 294(L) x 205(H) mm

Weight: 12 kg


Battery System

MESA-50 has its internal battery, whose life is up to 6 hours.



Option: Lightweight carrying case

Automatic measurement time extension




This function automatically sets up the optimum measurement time, based on judgment of results according to administrative values for the harmful element concentrations.

- #1 Measurement time varies with concentration of harmful element and material type of sample
- #2 The measurement will be terminated when OK or NG can be determined
- #3 If NG, set extended measurement time and repeat number (1-9 times)
- #4 Chemical analysis such as ICP would be recommended if it's impossible to determine OK/NG after maximum repeat measurement is finished at #2 .

MESA-50 specification

Basic items	Principle	Energy dispersive X-ray fluorescence spectrometry
	Target application	RoHS, ELV, Halogen-Free
	Measurement Elements	13Al - 92U
	Sample types	Solid, Liquid, Powder
X-ray generator	X-ray tube	Max 50kV, 0.2mA
	X-ray analysis diameter	1.2mm, 3mm, 7mm (Automatic switching)
	X-ray primary filter	4 types (Automatic switching)
Detector	Type	SDD (Silicon Drift Detector)
	Signal processor	DPP (Digital Pulse Processor)
Sample chamber	Atmosphere	Air
	Sample imaging	CCD camera
Utility	Operation	PC (Windows®7)
	Power supply	AC adapter (100-240V, 50/60Hz), battery
Software	Analysis Function	Multilayer Film FPM (Optional), Sb/As analysis (Optional)

Dimensional outline



Installation example



Options & Consumables



Lightweight carrying case



Cable fixing jig



Disposable cell



Transmission X-ray film



Cell window

* Windows and Excel are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries.



POWER

CHARGE

ERROR



Dedicated analyzers for WEEE/RoHS

MESA-50K

X-RAY FLUORESCENCE ANALYZER



Chamber size:
460(W) × 360(D) × 150(H) mm

New!
Big Chamber



XGT-5200WR

X-RAY ANALYTICAL MICROSCOPE



Chamber size:
350(W) × 400(D) × 40(H) mm



Please read the operation manual before using this product to assure safe and proper handling of the product.

- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.

<http://www.horiba.com> e-mail: info@horiba.co.jp

● **HORIBA, Ltd.**
Head Office
2 Miyanoigashi, Kisshoin
Minami-ku, Kyoto, Japan
Phone: 81 (75) 313-8123
Fax: 81 (75) 321-5725

● **HORIBA (China) Trading Co., Ltd.**
Shanghai Office
Unit D, 1F, Building A, Synnex
International Park, 1068
West Tianshan Road,
Shanghai, 200335 China
Phone: 86 (21) 6289-6060
Fax: 86 (21) 6289-5553

Beijing Branch
12F, Metropolis Tower, No.2,
Haidian Dong 3 Street,
Beijing, 100080, China
Phone: 86 (10) 8567-9966
Fax: 86 (10) 8567-9066

Guangzhou Branch
Room 1611/1612, Goldlion
Digital Network Center,
138 Tiyu Road East,
Guangzhou, 510620, P.R.C.
Phone: 86 (20) 3878-1883
Fax: 86 (10) 3878-1810

● **HORIBA Korea Ltd.**
10, Dogok-Ro, 6-Gil,
Gangnam-Gu, Seoul,
135-860, Korea
Phone: 82 (2) 753-7911
Fax: 82 (2) 756-4972

● **HORIBA India Private Limited**
246, Okhla Industrial Estate,
Phase 3 New Delhi - 110020, India
Phone: 91 (11) 4669 5001
Fax: 91 (11) 4669 5010

● **HORIBA Instruments (Singapore) Pte Ltd.**
Head Office
10, Ubi Crescent #05-12
Lobby B Ubi Techpark
Singapore 408564
Phone: 65 (6) 745-8300
Fax: 65 (6) 745-8155

● **PT HORIBA Indonesia**
Jl. Jalur Sutera Blok 20A,
No. 16-17, Kel. Kunciran,
Kecamatan Pinang
Tangerang-15144,
Indonesia
Phone: 62 (21) 3044 8525
Fax: 62 (21) 3044 8521

● **HORIBA (Thailand) Ltd.**
393, 395, 397, 399, 401,
403, Latya Road,
Somdetchaopraya,
Klongsan, Bangkok 10600,
Thailand
Phone: 66 (0) 2 861 5995
ext.123
Fax: 66 (0) 2 861 5200

● **HORIBA Instruments Inc.**
Edison Office
3880 Park Avenue, Edison,
NJ 08820, U.S.A.
Phone: 1 (732) 494-8660
Fax: 1 (732) 549-5125

Irvine Office
9755 Research Drive,
Irvine, CA 92618, U.S.A.
Phone: 1 (949) 250-4811
Fax: 1 (949) 250-0924

● **HORIBA JOBIN YVON SAS**
16-18 rue du Canal,
91165, Longjumeau Cedex,
France
Phone: 33 (1) 69-54-13-00
Fax: 33 (1) 69-09-07-21

Spain Branch
Avenida Caidos de la Division
Azul 16, 28016 Madrid, Spain
Phone: 84 (4) 3795-8552
Fax: 84 (4) 3795-8553

● **HORIBA JOBIN YVON GmbH**
Haupt-Str. 1* D-82008
Unterhaching, Germany
Phone: 49 (89) 46-23-17-0
Fax: 49 (89) 46-23-17-99

● **HORIBA JOBIN YVON Srl**
Via Luca Gaurico
209 - 00143, Roma, Italy
Phone: 39 (6) 51-59-22-1
Fax: 39 (6) 51-96-43-34

● **HORIBA UK Ltd.**
2 Dalston Gardens,
Stanmore, Middx HA7 1BQ,
Great Britain, UK
Phone: 44 (208) 204-8142
Fax: 44 (208) 204-6142

● **HORIBA Instruments Brazil, Ltda.**
Rua: Presbitero Plinio Alves de
Souza, 645, Loteamento Polo
Multivas Barirro Medeiros -
Jundiai Sao Paulo -
CEP 13.212-181 Brazil
Phone: 55 (11) 2923-5400
Fax: 55 (11) 2923-5490