






## Fluorolog<sup>®</sup> Series


Modular Research Fluorometers for Lifetime and Steady State Measurements



Most sensitive fluorometer.  
CCD options for fast EEMs and mapping

# Fluorolog<sup>®</sup> Series

*"The modular Fluorolog can grow with your needs, adding capabilities as your changing research calls for it."*




The Fluorolog®-3 series encourages you to custom-tailor a spectrofluorometer's performance to the work you need to accomplish. HORIBA's all-reflective optics and high powered Xenon arc lamp offer the highest sensitivity on the market, across all wavelengths. Adding TCSPC with the DeltaTime, brings the power of fluorescence and phosphorescence lifetimes to your instrument, with the fastest collection electronics available.

Fully modular and customizable, the Fluorolog-3 can stand alone as the centerpiece of your laboratory, or interface with other instruments to multiply the capabilities of your lab.

The latest generation of Fluorolog-3 can include optional T-side optics to collect fast anisotropy measurements, or simply more detectors. Our Nanolog™ configuration offers high speed multi-channel detection from UV to NIR and dedicated software for today's nanomaterials.

Fluorolog-3-based solutions are available for measuring solid and liquid samples, with high throughput screening, fast 3D scans, cryogenic or elevated temperatures, absolute quantum yields, microliter volumes, stopped-flow mixing or titration, and even micron scale measurements using a microscope. With the industry's most extensive list of accessories, the Fluorolog-3 series offers unparalleled flexibility to meet all of your lab's experimental needs.




# Modular Design

HORIBA's Fluorolog-3 is manufactured with our own gratings, CCDs, InGaAs array detectors, and fast PMTs, offering you everything you need to build the perfect system. Our quality components are housed in solidly constructed compartments with all-reflective optics to provide the highest light throughput and sensitivity on the market. Below are four classic examples of how to configure a system. DeltaTime can easily be added to any system to bring the fastest TCSPC collection to your laboratory.

## Fluorolog 3-11

**The basic configuration** consists of single-grating monochromators in excitation and emission positions, and a red-sensitive photomultiplier. Add any accessory now, or expand your capabilities later.

The FL-3-11 provides outstanding sensitivity and performance for most experimental requirements.




## Fluorolog 3-21

**One of our most popular configurations**, the Fluorolog 3-21 configuration is an L-Format fluorimeter with a double excitation monochromator on the excitation channel and a single emission monochromator on the right hand side of the sample compartment. This configuration offers excellent stray light rejection and versatility for most UV-Vis applications, and is ideal for highly scattering samples.



## Fluorolog 3-221

**The ultimate in stray-light rejection**, the double-grating monochromators in excitation and emission positions are perfect for highly scattering biological samples like lipids and proteins, or solids like powders, semiconductors, or phosphors. You also get a bonus in sensitivity. The additive grating design allows you to open the slits twice as wide for the same resolution you would get in a single-grating monochromator. Standard intermediate slits let you reduce stray-light even further. This system also adds a second emission monochromator, so-called T-format, for simultaneous dual emission detection.



## Nanolog

**Dedicated SWCNT and nanomaterials fluorometer.** Alternate between the best in scanning resolution and stray light rejection to the instantaneous acquisition and spatial resolution of an imaging spectrometer with a NIR InGaAs array and UV-Vis CCD. Our NanoLog is the prime example of this configuration, specially optimized for analyzing single-walled carbon nanotubes (SWCNT), quantum dots, and other nanomaterials.



# A look under the hood

Powerful 450 W Xe arc lamp  
for high intensity, broad  
spectrum excitation.


Pulsed Xe flash lamp  
for phosphorescence  
measurements

All-reflective optics give the  
highest throughput across  
all wavelengths

Near IR PMT for steady-state  
or TCSPC measurements.  
Choose from a number of PMTs,  
photodiodes, or array detectors to  
optimize your wavelength range,  
sensitivity, or collection speed.

T-side detection optics for  
simultaneous detection anisotropy  
or a second detection range.


Removable sample drawer to quickly exchange  
with thermostated holders, fiber launchers,  
multi-sample changers, or TCSPC light sources.



Double monochromators provide the highest stray light rejection, highest spectral resolution, and highest light throughput. Multi-channel spectrographs are available options.

Front face optics for efficient collection of light from solid samples, powders, films, and opaque solutions.

The industry standard R928P PMT. HORIBA gets specially-selected PMTs for our instruments to ensure the best performance.



# Considering Options

## **Single monochromator vs. double monochromator**

A double monochromator is recommended for solid (powders, films, devices, etc.) and turbid (cell suspensions, colloids, aggregates, etc.) samples, where scattered light can significantly interfere with the signal.

“Stray light” is unintended light from the sample or optical path that reaches the detector. It contributes to the background signal, and in extreme cases overwhelms the desired signal. The best strategy for removing stray light is the use of a second monochromator in series. HORIBA uses a true double monochromator with an adjustable center slit for the best stray light rejection available.

## **L- vs. T-format**

T-format is best for multiple detectors and real-time anisotropy usage.


The Fluorolog is modular and accessible, allowing the user to trade one detector for another. If your research will routinely switch between multiple detectors, a great deal of time will be saved with a dedicated second emission monochromator for your other detector. This is particularly important if you want to optimize one detection channel for Vis/UV, and the other for NIR.

In addition, if you plan to do many anisotropy measurements or anisotropy kinetics, T-format allows you to efficiently collect two polarization angles simultaneously. This decreases measurement time, and increases kinetics time resolution.

## **Array vs. PMT**

Array detectors are faster and will give greater S/N for a given collection time. PMTs are more sensitive for extremely weak samples, and necessary for any time-domain measurements

CCDs or InGaAs arrays collect a wide spectral range in a fraction of a second. If you have many samples to collect, want to collect excitation-emission matrices (EEMs), want to minimize photobleaching/photo-damage, or simply want to save time, you should consider an array detector. If, on the other hand, you expect to do phosphorescence or TCSPC measurements, a PMT is essential. The very weakest of samples will have a slight advantage with a PMT, and if this is the focus of your research, a PMT is your better choice.




### **PMT vs. Photodiodes**

Photomultiplier tubes provide better sensitivity and better repeatability. Photodiodes can cover very large ranges well into the NIR, and can be more affordable.


We use our PMTs in digital, single-photon counting mode, which means each collected photon is registered as a single count. The noise from this kind of measurement is simple and easily predicted, making single-photon counting the high-mark of fluorescence measurement for many years. In the near IR, we offer NIR PMT detectors and analog photodiodes. If you need detection beyond 1700 nm, photodiodes are the only option, and we can support measurements extending beyond 20  $\mu\text{m}$ .

### **Steady State vs. Lifetime**



Frequently, the information of a steady-state spectrum is sufficient, but the additional information of time-domain spectroscopy can clarify ambiguous data. Most fluorescence spectra are broad and often overlap each other, making quantification difficult. In biomolecules, there may be several Tryptophans fluorescing from different regions of a protein with very similar spectral properties. In these cases, it is helpful to look in the time domain, as fluorophores with similar spectral properties may have drastically different lifetimes. Furthermore, lifetime spectroscopy avoids photobleaching artifacts that confound steady-state measurements.



# Detector Options




# Optimize for the Near IR



• **Upconverters excited with HORIBA's 980 laser, miniature integrating sphere.** Upconverters can be excited with HORIBA's 980 nm laser, or your own custom source. Emission spectra (left) were collected from an integrating sphere accessory (right).

- **Single channel photodiodes**
- **Options from 0.8 – 40  $\mu\text{m}$**

HORIBA offers a number of photodiode detectors that cover the NIR-IR from 800 nm to 40  $\mu\text{m}$ , allowing you to use our industry leading instrument on modern semiconductors and nanomaterials.



**Range**

- **Steady state or TCSPC PMT**
- **950-1700 nm**
- **LN or TE cooling for ultra low dark counts**

We use selected TE-cooled Hamamatsu® near IR PMTs to offer the sensitivity of single photon counting for weakly emitting samples, as well as the versatility of ultrafast lifetime measurements.




**Versatility**

- **1024 channel InGaAs Array**
- **800-1700 nm (optional 800-2200 nm)**
- **LN or TE cooling for ultra low dark counts**


HORIBA makes its own TE- and LN-cooled CCDs and InGaAs array detectors. Two detectors can cover from 200 – 2300 nm, and collect complete spectra in a millisecond. High-throughput or excitation-emission matrix (EEM) measurements should utilize this technology.



**Speed**



The Raman band of ultra-pure water was measured in a standard 1 cm quartz cuvette. HORIBA qualifies its instruments with two measures of signal-to-noise: Root Mean Square (RMS) and First Standard Deviation (FSD) of light from its standard Xenon lamp source.



50 fM fluorescein in 0.01 N NaOH. Subtracting the dark counts and the signal from the solvent alone reveals the clean fluorescein spectrum at ultra-low concentration.

## Maximizing signal, Minimizing noise

It can be confusing to try and compare the sensitivity of one instrument to another with the multiple specifications and differing definitions that exist.

Sensitivity of a spectrofluorometer is generally reported as a minimum measurable concentration of a fluorophore (commonly Fluorescein), or the signal-to-noise ratio (S/N) of the Raman scattering peak of water. Most companies opt for the second, as it is a more readily verified number, especially in the field. Unfortunately, however, this specification can be confusing since two calculation methods prevail in the industry: First Standard Deviation (FSD) and Root Mean Square (RMS).

To avoid confusion, HORIBA Scientific quotes both specifications. And no matter how you look at it, the Fluorolog series has always offered the most sensitive modular research fluorometers available<sup>1</sup>. This is the result of increased signal due to the optimized design of the all reflective optics, the photon counting signal collection method, and the reduction of stray light by the choice of gratings and baffling.

## Measure smaller samples, detect smaller changes

### Big discoveries come from small changes.

Don't let your research get lost in the noise. Using the most sensitive fluorometer means having confidence in small changes in your data. These are the same small changes that drive curiosity and promote new discoveries.

### Dilute and small volume samples

Higher sensitivity means you can also measure more dilute samples, or need less sample to begin with. Given the cost and effort often required for sample preparation, this means you save both time and money.

## Sensitivity means speed


The real limit on how rapidly you can acquire data is not how fast you can scan the grating, but how long you have to measure at each wavelength to get an acceptable signal to noise. The higher sensitivity of the Fluorolog series means that not only can you measure weaker samples, but also that you can actually acquire good quality spectra faster.

1. Only the QuantaMaster 8000 series fluorometers, also from HORIBA Scientific, match the sensitivity of the Fluorolog.

## Dynamic Range

Measure more than 6 orders of magnitude of signal levels in the same spectrum


Fortunately not all samples are weak. In fact, some samples are quite strong. But some are both, which creates a challenge to measure the strong and weak peaks in the same scan. With an intra-spectral range of over 6 orders of magnitude, even wildly varying multi-component spectra are no problem to measure in one scan. Not only does this save time, but it is essential for kinetics studies when one cannot afford to repeat a measurement at different integration times to keep all peaks within range.



Spectrum of a mixture of Terbium Chloride and Europium Chloride in ultrapure water. The inset is a zoomed view of the spectrum showing the very weak Europium Chloride peak measured in the same scan as the extremely strong Terbium Chloride peaks recorded in the same measurement. The spectra of pure Terbium Chloride and pure Erbium Chloride are also shown in the zoomed in view for reference.

## Front Face

Solid and highly scattering liquid samples are more efficiently collected at front face. This gives maximal illumination of the sample. Double monochromators offer the best stray light rejection for highly scattering samples. The Fluorolog offers the best sensitivity not only due to high throughput of light, but due to its ability to reject stray light.



The plot of tryptophan in hemoglobin compares the spectra for right-angle and front-face viewing of a hemoglobin sample. For many years, hemoglobin was thought to be non-fluorescent because the fluorescence could not be detected at the conventional 90° angle. With front-face viewing, the fluorescence spectrum for the  $\beta 37$  tryptophan in the  $\alpha\beta$  region of hemoglobin can be easily obtained.

# CCD and InGaAs Array Detectors

HORIBA can offer the speed and performance of high-end array detectors at an affordable price, because we design and manufacture our own! These multi-channel detectors are ideal for high throughput, EEMs and microscopy mapping.

## The Multiplex Advantage:

### Low noise:

Our CCD's and InGaAs arrays are thermoelectrically or liquid nitrogen cooled, virtually eliminating their thermal noise. This helps the signal from your weakest samples stand out even more.


### Speed of acquisition:

Array detectors let you collect the entire spectrum in a fraction of a second. This facilitates possibilities of fast Excitation Emission Matrices (EEMs), full spectrum kinetics, or high throughput measurements of your samples.

### Signal to noise:


Simultaneous acquisition of the entire luminescence spectrum means more signal can be collected for a given period of time.

In a side-by-side comparison we took a 45 second measurement of Rhodamine B with both detectors under equal conditions. Even with 1/2 the wavelength resolution (PMT: 1 nm steps, CCD: 0.5 nm steps), the CCD shows 5x to 15x greater signal to noise. This is due to the better quantum efficiency of CCDs and the number of measurements that can be made with them in a short time.




50 repeated, 45 s emission scans of Rhodamine B were subtracted from the averaged scan to measure deviations. Measurements were adjusted for sample photobleaching. In the time given, the CCD (top) shows 5 to 15 times the signal to noise of the PMT (bottom).

# Absolute Quantum Yield Measurements




High precision absolute quantum yield characterization of a quantum dot sample enabled by the Quanta-Φ™ integrating sphere accessory.

Measure absolute quantum yields (QY) with the only accessory optimized for liquids, powders, films, and solids. The Quanta-Φ is based on a large 15 cm integrating sphere made of spectralon for the highest reflectivity over the broadest spectral range for the most accurate and reproducible QY values. The solid sample port can also accept wired sources for electroluminescence measurements.



FluoroLog-3 software also automatically generates CIE 1931 and 1976 values for your sample.

# Polarization




Fluorescein dissolved in four aqueous solutions of sucrose. With increased solution temperature, viscosity decreases, yielding faster rotation times and correspondingly lower anisotropies. Similarly, anisotropy is an excellent tool for understanding changes in macromolecule shape as well as molecular binding.

Adding optional polarizers to a FluoroLog-3 allows measuring anisotropy, and hence, the changes in molecular rotation in the sample. This indirect measure of the local viscosity gives information on sample aggregation, structural changes, molecular binding, and other mechanisms.

Flexibility

Enhanced Measurements

# Microvolumes



Using the Microsense, 5 µL of AlexaFluor 488 labeled IgG is enough to get steady-state or TCSPC data.

## Measure ultra-low sample volumes with easy sample recovery

HORIBA's Microsense, the most advanced microliter accessory on the market, is designed to get cuvette-quality steady state or lifetime data from a 1-5 µL of sample.


Avoiding dilution and nearly total recovery minimizes the need for sample, while preserving measurement sensitivity. Based on all quartz optics, Microsense is compatible with UV to NIR measurements.



# Well Plate Accessory

## Maximize your throughput

Search for "hits" or automate measurements of a large number of samples.



The MicroMax 384 accessory lets you automate your measurements for high-throughput data collection. Based on standard microtiter well plates, MicroMax saves you time and money by automating spectral, kinetics, single point, and time-resolved measurements on anywhere from 6 to 384 samples at a time. Fluorolog-3 software includes calibration curve routines for simple quantification of your results.

# Full Spectral Microscopy and FLIM


All instruments in the Fluorolog family can be fiber-coupled to virtually any upright or inverted microscope. Our confocal microscope coupling allows you to go beyond simple filter-based fluorescence microscopy, to full spectral analysis of spatially varying or extremely small volume samples.

Couple a microscope to a Fluorolog-3 and:

- Measure a complete spectrum of a sample as small as 1  $\mu\text{m}$ .
- Get spectra from as little as a few molecules of sample.
- Perform localized FRET measurements.


Add an automated stage and camera, and:

- Create complete spatial/spectral maps.
- Perform repetitive QC characterization of structured samples like photovoltaics.



Data from a real and forged US 5 dollar bill. Fluorescence imaging shows no obvious differences. However it is easy to spectrally distinguish the real from a fake bill.


## Absorbance/Transmittance Accessory



Reabsorption of fluorescence photons can occur even at moderate concentrations, and lead to increasing distortion of spectra with increasing concentration (left). The Fluorolog-3 software can correct this phenomenon, (known as the inner-filter effect), using the sample absorbance measured by the absorbance accessory (right).

Add absorbance to your fluorescence measurement.

Correct your fluorescence spectra for inner-filter effects in concentrated samples without dilution.



Flexibility

Enhanced Measurements

# Peltier-based Heated/Cooled Cuvette Holders


Rapid temperature control in single or multiple cuvette models with built-in magnetic stirring.

Precise temperature control for precise data on your protein folding, micellization, solubility, conformation, phase, and rotational transitions.


Rapidly vary sample temperature over a range of -40 to 110°C. Fluorolog-3 software also simplifies automating temperature dependence measurements, including complex ramps and profiles.

## Temperature Bath

An alternative to Peltier-based units. Better precision (0.01°C), range (-25 to 200°C) and long term stability, but with much slower temperature slewing.





Thermal unfolding of Bovine Serum Albumin in PBS and 1,8-ANS followed using three observables: the quenching of intrinsic Tryptophan fluorescence by water, the quenching of the 1,8-ANS, and the increased 2nd order scattered light which reports turbidity from protein aggregation.



# Economic LN Temperature Measurement Dewar Accessory

Cryogenic temperatures enable measurements of fine structure, enhanced phosphorescence, and rare conformations/states often not possible at room temperature. This low cost accessory readily permits measurements of samples at 77°K.



Temperature can add thermal broadening to fluorescence spectra and increase phosphorescence quenching. Spectra of Naphthalene dissolved in Methanol measured at room temperature (298°K) and in a liquid nitrogen dewar accessory (77°K). The low temperature spectrum reveals rich vibrational structure and longer wavelength phosphorescence. Phosphorescence peaks are also shown magnified 10x for clarity.

## Optional LN<sub>2</sub> and He Cryostats

For greater temperature range and sample type flexibility, the Fluorolog-3 supports various cryostats offering temperature control down to 4°K.




Flexibility

Temperature Control

## Other Accessories

The Fluorolog-3 series includes the most comprehensive line of accessories that enable researchers to extend the utility of their instrument to as many experiments as possible. The following is a partial list of accessories available, in addition to those previously discussed.

- Auto-titrator (injector) dual syringe, dual valve.
- Stopped flow rapid kinetics accessory.
- Solid sample holder designed for viewing front-face fluorescence of thin films, powders, pellets, paper, fibers, or microscopic slides. Variable alignment angle.
- 2-position sample holder with magnetic stirring bar.
- 4-position sample holder with magnetic stirring bar.
- 250  $\mu$ l reduced volume cell.
- 500  $\mu$ l cuvette 5x5 mm.
- 20  $\mu$ l HPLC flowcell.
- Fiber optic probe, bifurcated, randomized. Ideal for samples which cannot fit inside a standard sample chamber. Requires fiber adapter.
- External trigger accessory allows researchers to record time of sample modification during kinetic fluorescence experimentation. Also allows remote trigger activation to start experiment.
- Sealed water standard in scratch-proof housing for water Raman S/N verification.
- Emission correction factor kit.
- Excitation correction factor kit.
- Purge port, quartz windows for sample compartment for use with nitrogen purging.



# Specifications: Fluorolog-3

Optics	All reflective optics for high sensitivity at all wavelengths and for microsamples
Source	450 W CW Ozone-free xenon arc lamp (250-2500 nm)
Monochromators	Czerny-Turner design with plane gratings for optimized focus at all wavelengths and minimum stray light
Excitation grating	1200 groove/mm blazed at 330 nm
Emission grating	1200 groove/mm blazed at 500 nm
Bandpass	0 to 30 nm (single mono, 1200 gr/mm grating), continuously adjustable 0 to 15 nm (double mono, 1200 gr/mm grating), continuously adjustable
Wavelength Accuracy	± 0.5 nm
Integration Time	1ms – 160 s
Base detector	Photomultiplier R928P, spectral coverage 200 nm – 870 nm,
Reference Detector	UV enhanced silicon photodiode
Water Raman S/N	>30,000:1 RMS (15,000:1 FSD)
Dimensions (FL3-11)	77.8 cm (w) x 34.9 cm (h) x 102.2 cm (d)*
<b>Lifetime options:</b>	
TCSPC:	
Lifetime range with standard detector	<150 ps to >10 s
Lifetime range with optional PPD detector	<25 ps to >10 s
Lifetime range with suitable laser and detector	<5 ps to >10 s

\*Dimensions vary for each FL-3 configuration.

Flexibility

Accessories


# Take your measurements to the next dimension: Time

## Lifetime measurements:



# TCSPC

(Time Correlated Single Photon Counting)


- Robust data, independent of concentration and photobleaching.
- Lifetimes from 5 ps to seconds.
- Unique SpectraLEDs for highest efficiency phosphorescence measurements.
- TCSPC lifetimes, anisotropy, TRES, and kinetics.
- Working in the time domain removes the confounding influences of concentration and photobleaching.
- Add lifetime measurement to any Fluorolog with the DeltaTime™ TCSPC accessory. DeltaTime has the same footprint as a mouse pad, but is powerful enough to deliver 12 decades of lifetimes. With its industry-leading true 100 MHz system operation, DeltaTime offers TCSPC acquisition rates, with all decays being acquired in mere milliseconds, allowing for TCSPC lifetime kinetics of fast reactions. Its crystal-locked timing circuits never require recalibration. Select from our current catalog of over 70 compact pulsed light sources, with more being added all the time. And once you get your data, powerful DAS6 analysis software lets you choose among 9 fitting models.




Lifetimes of Rhodamine 6G in Methanol measured using optional ultrafast PPD detector. At high concentrations, self-quenching results from homodimers and trimers formation. Lifetimes as short as 75 ps are seen, as well as homo-FRET at lower concentrations.



Non-extensive distribution fits: 1,8-ANS exists in free solution and partitions into several disparate environments in Bovine Serum Albumin. Within each environment, a distribution of states exists with a corresponding distribution of lifetimes. FluoroLog-3 software not only offers standard discrete exponential fitting but also several energy transfer and distribution models, including the proprietary non-extensive distribution shown here.




Sample photobleaching corrupts kinetic data by adding an exponential term. Lifetimes are more robust than steady-state intensities. To support rapid kinetics, the FluoroLog-3 is capable of measuring a complete lifetime decay in as little as every 10 ms. Example shown is the binding of Curcumin to Serum Albumin.



## Phosphorescence

Measure phosphorescence spectra, lifetimes from microseconds to seconds. The Fluorolog-3 can be upgraded with an optional pulsed xenon flash lamp, enabling lifetime measurements down to 10  $\mu\text{s}$  with no additional electronics or detectors. Ideal for measuring lanthanide tagged samples or rare earth phosphors used in lighting applications.



Complex solutions like this mixture of Bovine Serum Albumin (BSA) and Terbium Chloride ( $\text{Tb}^{3+}$ ) can be challenging to interpret. Using the pulsed light source of the Phosphorescence option allows you to temporally “gate out” the BSA fluorescence, leaving only the  $\text{Tb}^{3+}$  phosphorescence. Inset shows the single wavelength phosphorescence decay of  $\text{Tb}^{3+}$  in this mixture; colored circles correspond to the gated spectra in the main figure.

# FluorEssence™ Software

**Simple enough for the occasional user; Powerful enough for the most elaborate experiments.**

## Fluorescence software that works like you do

- Efficiently develop your experimental method and then save it for future use.
- Data collection, analysis and report generation are easily streamlined.
- Full software control of accessories.
- Automate repetitive experiments with a built-in batch mode.
- All instrument calibration parameters are automatically applied per method.
- Of course, all experimental parameters are always saved along with the data file for comparison with previously collected data.


## Convert data to answers

Powerful processing and data-management tools of OriginPro™ include a complete suite of data reduction tools.

## A complete library of video tutorials to get you started

### Features

- Data views in workbook formats, keeping graphs, tables and notes together for each experiment
- Integrate, differentiate, or fit fluorescence data to Gaussian, Lorentzian, and custom curves
- Zooming and scaling
- Contour maps and profiles from 3D plots
- Peak finding
- Standard arithmetic
- 3D perspective
- Smoothing
- Deconvolution
- Excitation/emission correction
- Interpolation and extrapolation
- Blank subtraction
- Normalization
- PLQY calculator wizard (for use with Quanta- $\phi$ ™ accessory)




# Nanosizer Software


Single-walled carbon nanotubes come in a large variety of structures and diameters. With a sample of these nanotubes, the NanoLog can run an excitation-emission matrix scan. With the finished scan, you transfer these parameters to our exclusive Nanosizer software, which automatically determines the chirality and diameter from the spectral landscape. Below is an actual scan and complete analysis run on a NanoLog:

- 3-D spectral surface simulation
- Simultaneous analytical simulation of spectral surfaces
- Rapid preliminary scanning to recognize peaks and their shapes for easy model-fitting
- Complete, easy-to-edit model-parameter table for nanotube mixtures
- Nanotube species-recognition with editable library
- Nanotube species-recognition with user's analytical simulations
- Complete reports and charts in common spreadsheet format
- Optional "enhanced" fitting-engines for statistically robust simulations


A



B




Emission-excitation scan of a mixture of single-walled carbon nanotubes. Plot A is the spectral landscape recorded by the NanoLog; plot B is the structural assignment determined by the Nanosizer software. Diameters and colors of circles in plot B are related to peak-intensities in plot A.



Screen-capture from FluorEssence data acquisition and analysis software showing selected excitation and emission spectra centered on the intense peak emitted from the (8,6) species of carbon nanotube determined in Plots A and B.

# HORIBA FLUORESCENCE

Scientific



Our Roots Grow Deep

[www.horiba.com/fluorescence](http://www.horiba.com/fluorescence)  
[info.sci@horiba.com](mailto:info.sci@horiba.com)

**HORIBA**  
 Scientific

- USA:** HORIBA Instruments, Inc., 3880 Park Avenue, Edison, NJ 08820-3012 - Toll-free: +1-866-562-4698 - Tel: +1 732 494 8660 - Fax: +1 732 549 5125 - Email: [info.sci@horiba.com](mailto:info.sci@horiba.com)
- France:** HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel: +33 (0)1 69 74 72 00 - Fax: +33 (0)1 69 09 07 21 - Email: [info-sci.fr@horiba.com](mailto:info-sci.fr@horiba.com)
- Japan:** HORIBA Ltd., Tokyo Branch Office, 2-6, KandaAwaji-cho, Chiyoda-ku, Tokyo 101-0063, Japan - Tel: +81-(0)3 6206 4721 - Fax: +81 (0)3 6206 4730 - Email: [info-sci.jp@horiba.com](mailto:info-sci.jp@horiba.com)
- Germany:** HORIBA Jobin Yvon GmbH, Hauptstrasse 1, 82008 Unterhaching - Tel: +49 (0)89 4623 17-0 - Fax: +49 (0)89 4623 17-99 - Email: [info-sci.de@horiba.com](mailto:info-sci.de@horiba.com)
- Italy:** HORIBA Jobin Yvon Srl., Via Cesare Pavese 21, 20090 Opera (Milano) - Tel: +39 2 5760 3050 - Fax: +39 2 5760 0876 - Email: [info-sci.it@horiba.com](mailto:info-sci.it@horiba.com)
- UK:** HORIBA UK Ltd., 2 Dalston Gardens, Stanmore, Middlesex HA7 1BQ - Tel: +44 (0)20 8204 8142 - Fax: +44 (0)20 8204 6142 - Email: [info-sci.uk@horiba.com](mailto:info-sci.uk@horiba.com)
- China:** HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel: +86 (0)21 6289 6060 - Fax: +86 (0)21 6289 5553 - Email: [info-sci.cn@horiba.com](mailto:info-sci.cn@horiba.com)
- Brazil:** HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multivias, Bairro Medeiros, Jundiaí / SP, CEP 13.212-181 - Tel: +55 (0)11 2923 5400 - Fax: +55 (0)11 2923 5490 - Email: [infocientifica.br@horiba.com](mailto:infocientifica.br@horiba.com)
- Other:** Tel: +1 732 494 8660 - Email: [info.sci@horiba.com](mailto:info.sci@horiba.com)

HORIBA Scientific has a policy of continuous product development, and reserves the right to amend part numbers, descriptions and specifications without prior notice.