


## THEORETICAL EFFICIENCY CURVE

### Reference 530 34 ---

Grating Type	Plane Holographic Grating
Groove Density (gr/mm)	600
Spectral Range (nm)	700-1750
Blaze Wavelength (nm)	1000
Blaze Angle (deg)	
Configuration	



This efficiency curve should be used as reference only and does not imply any guarantee on actual performance. Actual performance of the grating will depend on several variables including the usage geometry and the coating. Efficiency curves labeled as "Theoretical" are calculated estimates of performance based on groove profile.