


Results for the **Aqualog** river water testing demonstration in England


This Video (www.youtube.com/watch?v=-fycBm9apSk) demonstrates the **Aqualog** in a real-life scenario, generating EEMs for CDOM (Coloured Dissolved Organic Matter) analysis in a fast and simple process.

This Video was shot on location at a river complex in South East England.




Upstream of the Sewage Works Outlet

Evidence of humic-like material.



In the vicinity of the Sewage Works Outlet

Evidence of both humic-like and protein-like material.



Downstream of the Sewage Works Outlet

Protein-like material weaker signal.

www.aqualog.com

HORIBA
Scientific

JOBIN YVON
Technology