


## SPECTRAL EFFICIENCY CURVE

### Reference 510 03 ---

Grating Type	Plane Mechanically Ruled Grating
Groove Density (gr/mm)	1800
Spectral Range (nm)	
Blaze Wavelength (nm)	630
Blaze Angle (deg)	34°32'
Configuration	Quasi-littrow

### 510 03 ---



This efficiency curve should be used as reference only and does not imply any guarantee on actual performance. Actual performance of the grating will depend on several variables including the usage geometry and the coating. Efficiency curves labeled as "Theoretical" are calculated estimates of performance based on groove profile.