

## PGM - PGS 200

### *Plane Grating Monograph 200*

### *Vacuum Monochromator and Spectrograph*



- ❑ 3 Interchangeable Aberration Corrected Gratings
- ❑ High Throughput
- ❑ Interchangeable Kinematic Grating Mounts
- ❑ Focal Length 250 mm
- ❑ Deviation 165°
- ❑ Spectral Wavelength Range 4-80 nm
- ❑ Spectral Resolution of few Angstrom
- ❑ Monochromator and Spectrograph
- ❑ MCP/CCD Detector adaptation
- ❑ 25 mm Flat Field
- ❑ Removable Entrance and Exit Slit
- ❑ High Vacuum Compatible
- ❑ Fixed or Adjustable Entrance and Exits Slits
- ❑ Stepping Motor Drive

**We are pleased to work with you to find a suitable solution to your requirements.**

#### **In North America:**

HORIBA Jobin Yvon Inc.  
3880 Park Avenue  
Edison, New Jersey  
08820, USA  
Tel: +1 732 494 8660 ext. 328  
Fax: +1 732 548 9725  
Email: [custom\\_gratings@jobinyvon.com](mailto:custom_gratings@jobinyvon.com)

#### **In France:**

HORIBA Jobin Yvon  
16-18, rue du Canal  
91165 Longjumeau cedex  
France  
Tel: +33 (0)1 64 54 13 00  
Fax: +33 (0)1 69 09 93 19  
Email: [info@jobinyvon.fr](mailto:info@jobinyvon.fr)

#### **Others areas:**

China, Beijing: +86 (0) 10 6849 2216  
China, Shanghai: +86 (0) 21 6447 9785  
Germany: +49 (0) 89 46 23 17 0  
Italy: +39 0 2 57 60 30 50  
Japan: +81 (0) 3 3861 8234  
U.K.: +44 (0) 20 8204 8142