

# Aqualog Version 3.6.11.20: Release Notes

---

## Product Release Information

Product	Aqualog Software
Release Number	3.6.11.20
Release Date	11/01/2016
Customer Support	For more information or support, please visit our website

## Introduction

This is a release to add features to Aqualog software version 3.6.11.20. Your feedback is important to us. Please help us make this the best product possible. Keep in mind that we are continuing to work on Aqualog and things may change in the future.

## What's New

1. Implemented an Add Contour Plots option to allow flexibility in data saved.
2. Added support and description of how to simulate a synchronous scan plot with the arbitrary line feature in the Origin Contour Profiles tool.
3. Installer includes updated users' manual in PDF format.

## Installation and Upgrade Notes

Patch can be applied to Aqualog 3.6.3 or newer versions.

## Recent Changes

Type	Tracking #	Description
feature	1	Implemented an Add Contour Plots option to allow flexibility in data saved.
feature	2	Support for 2D graph overlay tool.
Enhancement	3	Changes to make Sample Queue more robust
Enhancement	4	Changed maximum number of Sample Queue samples from 1000 to 5000.
Enhancement	5	Adjusted the size of the dialog containing Continue and Cancel button so that graphics are supported on systems using larger font
Enhancement	6	Changes for compatibility with Windows 10
Enhancement	7	Enhancements to stability of continuous CCD acquisitions.

Enhancement	8	Allows water Raman signal-to-noise ratio to be run multiple times for stability testing
Enhancement	9	Remembers last used folder when selecting/saving files.
feature	10	Added support of 0 nm excitation (0 order) for emission experiment to allow reflected measurements to be collected
Enhancement	11	Installer includes updated users' manual in PDF format.
Enhancement	12	Revisions to the Aqualog manual for instructions on lamp alignment and replacement.